```
In [ ]:
          cd /content/drive/My Drive/Xray report generation
         /content/drive/My Drive/Xray report generation
In [ ]:
        ls
         AttentionModel.ipynb
                                                          Medical.csv
         Basic_Model.ipynb
                                                          model_1_plot.png
         brucechou1983_CheXNet_Keras_0.3.0_weights.h5
                                                          NLMCXR png/
         Chexnetmodel/
                                                         'Old Files'/
         'Copy of AttentionModel.ipynb'
                                                          preprocessed.csv
         'Copy of SimpleModel.ipynb'
                                                          simple encoder decoder plot.png
         dup preprocessed.csv
                                                          SimpleModel.ipynb
         ecgen-radiology/
                                                          tokenizer1.pkl
         image feature vector1.pkl
                                                          Untitled0.ipynb
          indiana projections.csv
                                                          Untitled1.ipynb
         logs/
In [ ]:
        import re
         import pandas as pd
         import matplotlib.pyplot as plt
         import string
         from sklearn.utils import shuffle
         from sklearn.model selection import train test split
         from tqdm import tqdm
         import numpy as np
         import tensorflow as tf
         from nltk.translate.bleu_score import sentence_bleu
In [ ]: | df = pd.read_csv('Medical.csv')
```

In []: df.head()

Out[5]:		image_name	image_caption	comparison	indication	findings	impression
	0	CXR1_1_IM-0001- 3001.png,CXR1_1_IM- 0001-4001.png	xray chest pa and lateral	none	positive tb test	the cardiac silhouette and mediastinum size ar	normal chest x
	1	CXR10_IM-0002- 2001.png,CXR10_IM- 0002-1001.png	pa and lateral chest x	chest radiographs	male chest pain	the cardiomediastinal silhouette is within nor	no acute cardiopulmonary process
	2	CXR100_IM-0002- 1001.png,CXR100_IM- 0002-2001.png	chest v frontallateral pm	none	no indication	both lungs are clear and expanded heart and me	no active disease
	3	CXR1000_IM-0003- 1001.png,CXR1000_IM- 0003-2001	pa and lateral chest x	pa and lateral chest radiographs	male	there is increased opacity within the right up	increased opacity in the right upper lobe with
	4	CXR1001_IM-0004- 1001.png,CXR1001_IM- 0004-1002.png	chest v frontallateral pm	none	dyspnea subjective fevers arthritis immigrant 	interstitial markings are diffusely prominent	diffuse fibrosis no visible focal acute disease

```
In [ ]:
         df.shape
 Out[6]: (3851, 9)
In [ ]: df_projections = pd.read_csv('indiana_projections.csv')
          df_projections.head()
 In [ ]:
Out[11]:
              uid
                               filename
                                        projection
                1 1_IM-0001-4001.dcm.png
           0
                                            Frontal
           1
                1 1_IM-0001-3001.dcm.png
                                            Lateral
           2
                2 2_IM-0652-1001.dcm.png
                                            Frontal
                2 2_IM-0652-2001.dcm.png
                                            Lateral
                3 3_IM-1384-1001.dcm.png
                                            Frontal
 In [ ]:
          img_path = 'NLMCXR_png/'
```

```
In [ ]: #Creating structured data from raw xml files
        columns = ["image_1", "image_2", "impression"]
        df 1 = pd.DataFrame(columns = columns)
        columns = ["image 1", "image 2", "impression"]
        df dup = pd.DataFrame(columns = columns)
        no lateral = 0
        for item in tqdm(df.iterrows()):
            l = item[1]['image_name'].split(',')
            if len(1) > 2:
                 li, last_img = find_fr_lr_images(1)
                 if last img == "":
                     no lateral +=1
                     li, last_img = li[:-1], li[-1]
                for i in li:
                     image_1 = i
                     image 2 = last img
                     df 1 = df 1.append(pd.Series([image 1, image 2,item[1]['impression']]
            elif len(1) == 2:
                 image_1 = 1[0]
                 image 2 = 1[1]
                 df 1 = df 1.append(pd.Series([image 1, image 2, item[1]['impression']], i
            elif len(1) == 1:
                #creating duplicate dataframe separately to keep it in all dataset train
                df_dup = df_dup.append(pd.Series([1[0], 1[0],item[1]['impression']], inde
        print("Total Report without Lateral images {}".format(no_lateral))
```

1155it [00:03, 312.59it/s]/usr/local/lib/python3.6/dist-packages/ipykernel_laun cher.py:10: DeprecationWarning: The truth value of an empty array is ambiguous. Returning False, but in future this will result in an error. Use `array.size > 0` to check that an array is not empty.

Remove the CWD from sys.path while we load stuff.
3851it [00:13, 291.25it/s]

Total Report without Lateral images 1

```
In [ ]:
           df 1.shape
Out[15]: (3532, 3)
 In [ ]:
           df 1.head()
Out[16]:
                                image_1
                                                          image_2
                                                                                                 impression
                CXR1_1_IM-0001-3001.png
                                          CXR1_1_IM-0001-4001.png
                                                                                               normal chest x
            1
                 CXR10 IM-0002-2001.png
                                           CXR10 IM-0002-1001.png
                                                                              no acute cardiopulmonary process
                CXR100 IM-0002-1001.png
                                           CXR100 IM-0002-2001.png
                                                                                            no active disease
            2
                       CXR1000 IM-0003-
                                                 CXR1000_IM-0003-
            3
                                                                    increased opacity in the right upper lobe with...
                                                          2001.png
                                1001.png
                       CXR1000 IM-0003-
                                                 CXR1000 IM-0003-
            4
                                                                    increased opacity in the right upper lobe with...
                                3001.png
                                                          2001.png
 In [ ]:
           df dup.shape
Out[17]:
           (446, 3)
           df dup.head()
 In [ ]:
Out[18]:
                                image_1
                                                          image_2
                                                                                                 impression
                       CXR1003 IM-0005-
                                                 CXR1003 IM-0005-
                                                                              retrocardiac soft tissue density the
            0
                                2002.png
                                                          2002.png
                                                                                                appearanc...
                       CXR1012 IM-0013-
                                                 CXR1012 IM-0013-
            1
                                                                    bibasilar airspace disease and bilateral pleur...
                                1001.png
                                                          1001.png
                       CXR1024 IM-0019-
                                                 CXR1024 IM-0019-
            2
                                                                                         no acute abnormality
                                1001.png
                                                          1001.png
                       CXR1026 IM-0021-
                                                 CXR1026 IM-0021-
            3
                                                                              no acute cardiopulmonary disease
                                2002.png
                                                          2002.png
                       CXR1029 IM-0022-
                                                 CXR1029 IM-0022-
            4
                                                                        no pneumonia heart size normal scoliosis
                                1001.png
                                                          1001.png
           Adding Start and end tokens
 In [ ]:
          def add_start_end_token(data):
                # Combining all the above stundents
                preprocessed reviews eng = []
                # tqdm is for printing the status bar
                for sentance in tqdm(data.values):
                     sentance = '<start> ' + sentance + ' <end>'
                     preprocessed reviews eng.append(sentance.strip())
                return preprocessed reviews eng
```

```
In [ ]: | df 1['impression'] = add start end token(df 1['impression'].astype(str))
                         3532/3532 [00:00<00:00, 771697.75it/s]
 In [ ]:
           df dup['impression'] = add start end token(df dup['impression'].astype(str))
                              446/446 [00:00<00:00, 396242.23it/s]
 In [ ]:
           df 1.head()
Out[22]:
                                image_1
                                                          image_2
                                                                                                  impression
               CXR1 1 IM-0001-3001.png
                                          CXR1 1 IM-0001-4001.png
                                                                                   <start> normal chest x <end>
                                                                        <start> no acute cardiopulmonary process
                 CXR10_IM-0002-2001.png
                                           CXR10_IM-0002-1001.png
            1
                                                                                                       <end>
                CXR100 IM-0002-1001.png
            2
                                          CXR100 IM-0002-2001.png
                                                                                 <start> no active disease <end>
                      CXR1000 IM-0003-
                                                CXR1000 IM-0003-
            3
                                                                     <start> increased opacity in the right upper I...
                                1001.png
                                                          2001.png
                      CXR1000 IM-0003-
                                                CXR1000 IM-0003-
                                                                     <start> increased opacity in the right upper I...
                               3001.png
                                                          2001.png
           df dup.head()
 In [ ]:
Out[23]:
                                image_1
                                                          image_2
                                                                                                  impression
                      CXR1003 IM-0005-
                                                CXR1003 IM-0005-
            0
                                                                      <start> retrocardiac soft tissue density the a...
                               2002.png
                                                         2002.png
                      CXR1012 IM-0013-
                                                 CXR1012 IM-0013-
            1
                                                                     <start> bibasilar airspace disease and bilater...
                                1001.png
                                                          1001.png
                      CXR1024 IM-0019-
                                                CXR1024 IM-0019-
            2
                                                                              <start> no acute abnormality <end>
                                1001.png
                                                          1001.png
                      CXR1026 IM-0021-
                                                CXR1026 IM-0021-
                                                                        <start> no acute cardiopulmonary disease
            3
                               2002.png
                                                          2002.png
                                                                                                       <end>
                      CXR1029 IM-0022-
                                                 CXR1029 IM-0022-
                                                                          <start> no pneumonia heart size normal
            4
                                1001.png
                                                          1001.png
                                                                                                     scolios...
           df 1.to csv('preprocessed.csv')
           df dup.to csv('dup preprocessed.csv')
```

```
In [ ]: | df 1.info()
          <class 'pandas.core.frame.DataFrame'>
          RangeIndex: 3532 entries, 0 to 3531
          Data columns (total 3 columns):
           #
                Column
                              Non-Null Count
                                                Dtype
           - - -
           0
                              3532 non-null
                                                object
                image 1
                              3532 non-null
                                                object
           1
                image 2
                             3532 non-null
            2
                impression
                                                object
          dtypes: object(3)
          memory usage: 82.9+ KB
In [ ]:
          df dup.info()
          <class 'pandas.core.frame.DataFrame'>
          RangeIndex: 446 entries, 0 to 445
          Data columns (total 3 columns):
           #
                Column
                              Non-Null Count
                                                Dtype
           0
                              446 non-null
                                                object
                image 1
                              446 non-null
                                                object
           1
                image_2
            2
                impression 446 non-null
                                                object
          dtypes: object(3)
          memory usage: 10.6+ KB
          df 1 = pd.read csv('preprocessed.csv')
In [ ]:
          df 1.shape
In [ ]:
Out[9]: (3532, 4)
          df 1.head()
In [ ]:
Out[10]:
              Unnamed:
                                     image_1
                                                           image_2
                                                                                           impression
                                                   CXR1_1_IM-0001-
                             CXR1 1 IM-0001-
           0
                      0
                                                                              <start> normal chest x <end>
                                     3001.png
                                                          4001.png
                              CXR10 IM-0002-
                                                    CXR10_IM-0002-
                                                                          <start> no acute cardiopulmonary
           1
                      1
                                     2001.png
                                                          1001.png
                                                                                         process <end>
                                                   CXR100_IM-0002-
                             CXR100 IM-0002-
                      2
           2
                                                                            <start> no active disease <end>
                                     1001.png
                                                          2001.png
                            CXR1000 IM-0003-
                                                  CXR1000 IM-0003-
                                                                        <start> increased opacity in the right
           3
                      3
                                     1001.png
                                                          2001.png
                                                                                              upper I...
                                                  CXR1000 IM-0003-
                            CXR1000 IM-0003-
                                                                        <start> increased opacity in the right
                      4
                                                          2001.png
                                     3001.png
                                                                                              upper I...
 In [ ]: df 2 = pd.read csv("dup preprocessed.csv")
```

```
In [ ]: df_2.head()
```

```
Out[12]:
                Unnamed:
                                           image_1
                                                                    image_2
                                                                                                          impression
                         0
                                 CXR1003 IM-0005-
                                                          CXR1003 IM-0005-
                                                                               <start> retrocardiac soft tissue density the
             0
                         0
                                           2002.png
                                                                    2002.png
                                                                                  <start> bibasilar airspace disease and
                                 CXR1012 IM-0013-
                                                          CXR1012 IM-0013-
             1
                         1
                                                                    1001.png
                                           1001.png
                                                                                                              bilater...
                                 CXR1024 IM-0019-
                                                          CXR1024 IM-0019-
             2
                         2
                                                                                    <start> no acute abnormality <end>
                                           1001.png
                                                                    1001.png
                                 CXR1026 IM-0021-
                                                          CXR1026 IM-0021-
                                                                               <start> no acute cardiopulmonary disease
             3
                         3
                                                                    2002.png
                                           2002.png
                                 CXR1029 IM-0022-
                                                          CXR1029 IM-0022-
                                                                                <start> no pneumonia heart size normal
                                                                    1001.png
                                           1001.png
                                                                                                             scolios...
```

```
In [ ]:
         image names = []
         for img in tqdm(df['image name'].str.split(',')):
           for i in range(len(img)):
             image names.append(img[i])
                | 3851/3851 [00:00<00:00, 894713.60it/s]
         image_names[0:4]
In [ ]:
Out[14]: ['CXR1_1_IM-0001-3001.png',
           'CXR1_1_IM-0001-4001.png',
          'CXR10 IM-0002-2001.png',
          'CXR10 IM-0002-1001.png']
In [ ]:
         import tensorflow as tf
         from tensorflow.keras.applications import densenet
         from tensorflow.keras.applications.densenet import preprocess_input
         from tensorflow.keras.layers import Dense, Dropout,Input,Conv2D
         from tensorflow.keras.models import Model
In [ ]: chex = densenet.DenseNet121(include top=False,weights=None,input shape=(224,224,3
         X = chex.output
         X = Dense(14,activation='sigmoid',name='predictions')(X)
         model = Model(inputs=chex.input,outputs=X)
         model.load_weights('brucechou1983_CheXNet_Keras_0.3.0_weights.h5')
         image features model = Model(inputs=model.input,outputs= model.layers[-2].output)
 In [ ]:
         # image features model.summary()
```

```
In [ ]: | ## the input layer shape is (None, None, None, 3)
         ## the output layer is (None, 2048)
In [ ]: | #### This cell takes around 2 hrs to run.
         # image feature vectors = []
         # for img in tqdm(image names):
             image = tf.io.read file(img path+str(img))
             image = tf.image.decode jpeq(image,channels=3)
             image = tf.image.resize(image,(299,299))
             image = preprocess input(image)
             image features = image features model(tf.constant(image)[None,:]) # here we a
             image features = tf.reshape(image features,[-1,image features.shape[1]])
             image feature vectors.append(image features)
In [ ]: | # print(image feature vectors[1])
        # Len(image feature vectors)
In [ ]:
In [ ]: # import pickle
         # with open('image feature vector.pkl','wb') as f:
             pickle.dump(image feature vectors,f)
 In [ ]: |
         import pickle
         with open('image feature vector1.pkl','rb') as f:
           image feature vector 1 = pickle.load(f)
In [ ]: | print(image_feature_vector_1[1])
         tf.Tensor(
         [[5.0062913e-04 2.0340595e-03 1.2797897e-03 ... 9.1456938e-01
           9.2038012e-01 7.4285048e-01]], shape=(1, 1024), dtype=float32)
         Train, Validation and Test split
         i train,input test,o train,output test = train test split(df 1[['image 1','image
         input train,input validate,output train,output validate = train test split(i trai
In [ ]: input train.shape,input validate.shape,input test.shape,output train.shape,output
Out[20]: ((2542, 2), (636, 2), (354, 2), (2542,), (636,), (354,))
         i train dup,input test dup,o train dup,output test dup = train test split(df 2[['
         input_train_dup,input_validate_dup,output_train_dup,output_validate_dup = train_t
```

```
In [ ]: input train dup.shape,input_validate_dup.shape,input_test_dup.shape,output_train_
Out[22]: ((320, 2), (81, 2), (45, 2), (320,), (81,), (45,))
         Appending the duplicate data and the original data
In [ ]: | i train = np.append(input train,input train dup,axis=0)
         o_train = np.append(output_train,output_train_dup,axis=0)
         i validate = np.append(input validate,input validate dup,axis=0)
         o validate = np.append(output validate,output validate dup,axis=0)
         i test = np.append(input test,input test dup,axis=0)
         o test = np.append(output test,output test dup,axis=0)
In [ ]: | i_train[0]
Out[24]: array(['CXR914 IM-2417-1001.png', 'CXR914 IM-2417-3001.png'], dtype=object)
In [ ]: #https://datascience.stackexchange.com/questions/24511/why-should-the-data-be-shu
         # Here we will be shuffling the data
         for i in range(5):
           i_train,o_train = shuffle(i_train,o_train,random_state=15)
           i validate, o validate = shuffle(i validate, o validate, random state=15)
           i test, o test = shuffle(i test, o test, random state=15)
         TEXT TOKENIZATION
 In [ ]: from tensorflow.keras.preprocessing.text import Tokenizer
         from tensorflow.keras.preprocessing.sequence import pad sequences
         max_len_output = 60
         tokenizer = Tokenizer(oov token="\langle unk \rangle", filters='!"#$%&()*+.,-/:;=?@[\]^ `{|}~
         tokenizer.fit on texts(o train)
         text train = tokenizer.texts to sequences(o train)
         text test = tokenizer.texts to sequences(o test)
         text_val = tokenizer.texts_to_sequences(o_validate)
         dictionary = tokenizer.word index
         word2idx = \{\}
         idx2word = {}
         for k, v in dictionary.items():
             word2idx[k] = v
             idx2word[v] = k
In [ ]: vocab size = len(word2idx)+1
         vocab size
Out[27]: 1339
```

```
print("===== Top 6 Word and its Index =====")
         list(dictionary.items())[:6]
         ==== Top 6 Word and its Index =====
Out[28]: [('<unk>', 1),
          ('<start>', 2),
          ('<end>', 3),
          ('no', 4),
          ('acute', 5),
          ('cardiopulmonary', 6)]
In [ ]: ## Padding the text
         text train output = pad sequences(text train, maxlen=60, dtype='int32', padding='pos
         text validation output = pad sequences(text val,maxlen=60,dtype='int32',padding='
         text test output = pad sequences(text test, maxlen=60, dtype='int32', padding='post'
In [ ]: | text_train_output.shape
Out[30]: (2862, 60)
         Creating Tensorflow Dataset
 In [ ]: | # converting the numpy image files into tensorformats
         def multi image(img,imp):
           This function will be taking in the imgaes we have and converting into tensor for
            return tf.convert to tensor([image feature vector 1[image names.index(img[0].de
         dataset train = tf.data.Dataset.from tensor slices((i train,text train output))
In [ ]:
In [ ]:
         dataset validation = tf.data.Dataset.from tensor slices((i validate,text validation))
         dataset_train = dataset_train.map(lambda item1,item2: tf.numpy_function(multi_ima)
         dataset validation = dataset validation.map(lambda item1,item2: tf.numpy function
         BATCH SIZE = 32
In [ ]:
         BUFFER SIZE = 500
         embedding dimension = 256
         units = 512 # output
```

```
In [ ]: # input_shape is (2,1,1024)
# vocab_size = 1225
```

```
In []: vocab_size = 1339
    input_layer = tf.keras.layers.Input(shape=(2,1,1024))
    encoder_concat = tf.keras.layers.Concatenate()([input_layer[:,0], input_layer[:,1]
    encoder_out = tf.keras.layers.Dense(embedding_dimension)(encoder_concat)
    input_layer_text = tf.keras.layers.Input(shape=(1,))
    x = tf.keras.layers.Embedding(vocab_size, embedding_dimension)(input_layer_text)
    x = tf.keras.layers.Concatenate()([x, encoder_out])
    x = tf.keras.layers.LSTM(units)(x)
    x = tf.keras.layers.Dense(vocab_size)(x)
    model = tf.keras.Model([input_layer, input_layer_text],x)
    model.summary()
```

Model: "functional_5"

Layer (type)	Output Shape	Param # =======	Connected to
====== input_2 (InputLayer)	[(None, 2, 1, 1024)]	0	
 tf_op_layer_strided_slice (Tens	[(None, 1, 1024)]	0	input_2[0][0]
 tf_op_layer_strided_slice_1 (Te	[(None, 1, 1024)]	0	input_2[0][0]
input_3 (InputLayer)	[(None, 1)]	0	
concatenate (Concatenate) rided_slice[0][0]	(None, 1, 2048)	0	<pre>tf_op_layer_st tf_op_layer_st</pre>
rided_slice_1[0][0] 			
embedding (Embedding)	(None, 1, 256)	342784	input_3[0][0]
dense (Dense) [0]	(None, 1, 256)	524544	concatenate[0]
concatenate_1 (Concatenate)	(None, 1, 512)	0	embedding[0]
[0]			dense[0][0]
lstm (LSTM) [0][0]	(None, 512)	2099200	concatenate_1
dense_1 (Dense)	(None, 1339)	686907	lstm[0][0]
=======================================			

Total params: 3,653,435

Trainable params: 3,653,435 Non-trainable params: 0

```
In [ ]: tf.keras.utils.plot model(model, to file='simple encoder decoder plot.png', show sh
Out[40]:
                                                                                    input:
                                                                                            [(?, 2, 1, 1024)]
                                                                input_2: InputLayer
                                                                                    output: [(?, 2, 1, 1024)]
                                                                     (?, 2, 1, 1024)
                                                                                                                                             (?, 2, 1, 1024)
                                                             input:
                                                                                                                                     input:
                  tf_op_layer_strided_slice: TensorFlowOpLayer
                                                                                       tf\_op\_layer\_strided\_slice\_1: TensorFlowOpLayer
                                                             output:
                                                                      [(?, 1, 1024)]
                                                                                                                                     output:
                                                                                                                                              [(?, 1, 1024)]
                                                                                  input:
                                                                                          [(?, 1, 1024), (?, 1, 1024)]
                                                                                                                                           input:
                                                                                                                                                  [(?, 1)]
                                                          concatenate: Concatenate
                                                                                                                      input_3: InputLayer
                                                                                                (?, 1, 2048)
                                                                                  output:
                                                                                                                                          output: [(?, 1)]
                                                                                       input:
                                                                                               (?, 1, 2048)
                                                                                                                                                   (?, 1)
                                                                                                                 embedding: Embedding
                                                                         dense: Dense
                                                                                       output:
                                                                                                (?, 1, 256)
                                                                                                                                        output: (?, 1, 256)
                                                                                                                     [(?, 1, 256), (?, 1, 256)]
                                                                                                             input:
                                                                                  concatenate_1: Concatenate
                                                                                                             output:
                                                                                                                           (?, 1, 512)
                                                                                                            input:
                                                                                                                    (?, 1, 512)
                                                                                               lstm: LSTM
                                                                                                            output:
                                                                                                                     (?, 512)
                                                                                                              input:
                                                                                                                       (?, 512)
                                                                                              dense_1: Dense
                                                                                                              output: (?, 1339)
```

```
In [ ]: class Decoder(tf.keras.Model):
            def __init__(self, embedding_dim, units, vocab_size):
                 super(Decoder, self).__init__()
                 self.units = units
                 self.embedding = tf.keras.layers.Embedding(vocab size, embedding dim)
                 self.lstm = tf.keras.layers.LSTM(self.units,
                                                return sequences=True,
                                                return state=True,
                                                recurrent initializer=tf.keras.initializer
                self.dense = tf.keras.layers.Dense(vocab_size, kernel_initializer=tf.kera
            def call(self, x, features):
                #input x = input word teach forcing
                #input features = encoder image features
                x = self.embedding(x)
                x = tf.concat([x, tf.expand_dims(features,1)], axis=-1)
                output, state, _ = self.lstm(x)
                x = self.dense(output)
                 return x
```

Encoder Decoder Model

```
In [ ]: # modelling loss and accuracy for gradient tape based on this https://www.tensorf
        optimizer = tf.keras.optimizers.Adam()
        loss_obj = tf.keras.losses.SparseCategoricalCrossentropy(from_logits=True, reduc
                  = tf.keras.metrics.SparseCategoricalAccuracy()
        acc obj
        def loss_func(real,pred):
          # getting the mask
          mask = tf.math.logical not(tf.math.equal(real,0))
          # calculate the loss
          loss = loss obj(real,pred)
          # converting mask dtype to loss dtype
          mask = tf.cast(mask,dtype=loss .dtype)
          # applying the mask to loss
          loss = loss * mask
          loss = tf.reduce mean(loss )
          return loss
        def acc func(real, pred):
          acc f = acc obj(real,pred)
          return tf.reduce_mean(acc_f)
```

```
In [ ]: import datetime
    current_time = datetime.datetime.now().strftime("%Y%m%d-%H%M%S")
    train_log_dir = 'logs/Basic_gradient_tape/' + current_time +'/train'
    validation_log_dir = 'logs/Basic_gradient_tape/'+ current_time +'/test'
    train_summary_writer = tf.summary.create_file_writer(train_log_dir)
    validattion_summary_writer = tf.summary.create_file_writer(validation_log_dir)
```

```
In [ ]: encoder = Encoder(embedding_dimension)
  decoder = Decoder(embedding_dimension,units,vocab_size)
```

```
In [ ]:
        @tf.function
        def train step(tensor, target):
            loss = 0
            accuracy = 0
            dec input = tf.expand dims([tokenizer.word index['<start>']] * target.shape[0
            with tf.GradientTape() as tape:
                features = encoder(tensor)
                for i in range(1, target.shape[1]):
                    # passing the features through the decoder
                    predictions = decoder(dec_input, features)
                    loss += loss func(target[:, i], predictions)
                    accuracy += acc_func(target[:, i], predictions)
                    # using teacher forcing
                    dec input = tf.expand dims(target[:, i],1)
                    #print("decoder input teacher", dec_input.shape)
            total loss = (loss / int(target.shape[1]))
            total_acc = (accuracy / int(target.shape[1]))
            trainable variables = encoder.trainable variables + decoder.trainable variable
            gradients = tape.gradient(loss, trainable variables)
            optimizer.apply gradients(zip(gradients, trainable variables))
            return loss, total loss, total acc
        #validation function
        @tf.function
        def val step(tensor, target):
            loss val = 0
            accuracy_val = 0
            dec_input = tf.expand_dims([tokenizer.word_index['<start>']] * target.shape[0
            with tf.GradientTape() as tape:
                features = encoder(tensor)
                for i in range(1, target.shape[1]):
                    # passing the features through the decoder
                    predictions val = decoder(dec input, features)
                    loss_val += loss_func(target[:, i], predictions_val)
                    accuracy val += acc func(target[:, i], predictions val)
                    # using teacher forcing
                    dec_input = tf.expand_dims(target[:, i],1)
                    #print("decoder input teacher", dec_input)
            total loss val = (loss val / int(target.shape[1]))
            total acc val = (accuracy val / int(target.shape[1]))
            return loss val, total loss val, total acc val
```

MODEL Training

```
In [ ]: tf.keras.backend.clear session()
        EPOCHS = 10
        loss plot train = []
        loss plot val = []
        for epoch in range(0, EPOCHS):
            print("===== Start Epoch " +str(epoch + 1)+ " =======")
            total loss train = 0
            total acc train = 0
            total_loss_val = 0
            total acc val = 0
            print('Batchwise Train loss')
            for (batch, (jpg_tensor, target)) in enumerate(dataset_train):
                batch loss, t loss, t acc = train step(jpg tensor, target)
                total loss train += t loss
                total_acc_train += t_acc
                if batch % 40 == 0:
                    print ('Epoch {} Batch {} Loss {:.4f} acc {:.4f}'.format(
                       epoch + 1, batch, batch loss / int(target.shape[1]), t acc))
            loss_plot_train.append(total_loss_train / int(len(i_train) // BATCH_SIZE))
            with train summary writer.as default():
                tf.summary.scalar('loss', total_loss_train/ int(len(i_train) // BATCH_SIZ
                tf.summary.scalar('accuracy', total_acc_train/ int(len(i_train) // BATCH_
            print('Batchwise validation loss')
            for (batch, (jpg_tensor, target)) in enumerate(dataset_validation):
                batch loss val, t loss val, t acc val = val step(jpg tensor, target)
                total loss val += t loss val
                total_acc_val += t_acc_val
                if batch % 40 == 0:
                    print ('Epoch {} Batch {} Loss {:.4f} acc {:.4f}'.format(
                       epoch + 1, batch, batch_loss_val / int(target.shape[1]), t_acc_val)
            with validattion_summary_writer.as_default():
                tf.summary.scalar('loss', total loss val/int(len(i validate) // BATCH SIZ
                tf.summary.scalar('accuracy', total_acc_val/int(len(i_validate) // BATCH_
            template = 'Epoch {}, Loss: {}, Accuracy: {}, Test Loss: {}, Test Accuracy: {
            print (template.format(epoch+1,
                                     total loss train/ int(len(i train) // BATCH SIZE),
                                     (total_acc_train/ int(len(i_train) // BATCH_SIZE))*10
                                     total loss val/int(len(i validate) // BATCH SIZE),
                                     (total acc val/int(len(i validate) // BATCH SIZE))*10
```

```
41300559043884, Test Accuracy: 82.40718078613281
===== Start Epoch 2 ======
Batchwise Train loss
Epoch 2 Batch 0 Loss 0.5559 acc 0.7840
Epoch 2 Batch 40 Loss 0.4391 acc 0.7782
Epoch 2 Batch 80 Loss 0.3942 acc 0.7741
Batchwise validation loss
Epoch 2 Batch 0 Loss 0.4417 acc 0.7728
Epoch 2, Loss: 0.4169798791408539, Accuracy: 78.60909271240234, Test Loss: 0.39
965489506721497, Test Accuracy: 80.63554382324219
===== Start Epoch 3 ======
Batchwise Train loss
Epoch 3 Batch 0 Loss 0.3044 acc 0.7703
Epoch 3 Batch 40 Loss 0.2925 acc 0.7681
Epoch 3 Batch 80 Loss 0.2865 acc 0.7662
Batchwise validation loss
Epoch 3 Batch 0 Loss 0.3120 acc 0.7661
Epoch 3, Loss: 0.33620578050613403, Accuracy: 77.6750259399414, Test Loss: 0.35
80491840839386, Test Accuracy: 80.03739929199219
===== Start Epoch 4 ======
Batchwise Train loss
Epoch 4 Batch 0 Loss 0.2925 acc 0.7647
Epoch 4 Batch 40 Loss 0.3061 acc 0.7643
Epoch 4 Batch 80 Loss 0.3507 acc 0.7629
Batchwise validation loss
Epoch 4 Batch 0 Loss 0.2451 acc 0.7629
Epoch 4, Loss: 0.3041955232620239, Accuracy: 77.23392486572266, Test Loss: 0.33
000192046165466, Test Accuracy: 79.68699645996094
===== Start Epoch 5 ======
Batchwise Train loss
Epoch 5 Batch 0 Loss 0.2505 acc 0.7619
Epoch 5 Batch 40 Loss 0.3086 acc 0.7613
Epoch 5 Batch 80 Loss 0.2852 acc 0.7609
Batchwise validation loss
Epoch 5 Batch 0 Loss 0.2808 acc 0.7608
Epoch 5, Loss: 0.2833741307258606, Accuracy: 76.99518585205078, Test Loss: 0.31
85734152793884, Test Accuracy: 79.48360443115234
===== Start Epoch 6 ======
Batchwise Train loss
Epoch 6 Batch 0 Loss 0.3131 acc 0.7600
Epoch 6 Batch 40 Loss 0.3046 acc 0.7601
Epoch 6 Batch 80 Loss 0.1760 acc 0.7595
Batchwise validation loss
Epoch 6 Batch 0 Loss 0.3143 acc 0.7594
Epoch 6, Loss: 0.27113163471221924, Accuracy: 76.8564453125, Test Loss: 0.31030
49397468567, Test Accuracy: 79.36990356445312
===== Start Epoch 7 ======
Batchwise Train loss
Epoch 7 Batch 0 Loss 0.3416 acc 0.7588
Epoch 7 Batch 40 Loss 0.3698 acc 0.7588
Epoch 7 Batch 80 Loss 0.2490 acc 0.7584
Batchwise validation loss
Epoch 7 Batch 0 Loss 0.3573 acc 0.7584
Epoch 7, Loss: 0.26025623083114624, Accuracy: 76.72484588623047, Test Loss: 0.3
112887442111969, Test Accuracy: 79.25816345214844
===== Start Epoch 8 ======
Batchwise Train loss
```


```
Epoch 8 Batch 0 Loss 0.1634 acc 0.7581
        Epoch 8 Batch 40 Loss 0.2906 acc 0.7582
        Epoch 8 Batch 80 Loss 0.3026 acc 0.7582
        Batchwise validation loss
        Epoch 8 Batch 0 Loss 0.4076 acc 0.7578
        Epoch 8, Loss: 0.25657811760902405, Accuracy: 76.65787506103516, Test Loss: 0.3
        164512813091278, Test Accuracy: 79.20201110839844
        ===== Start Epoch 9 ======
        Batchwise Train loss
        Epoch 9 Batch 0 Loss 0.2662 acc 0.7575
        Epoch 9 Batch 40 Loss 0.2277 acc 0.7575
        Epoch 9 Batch 80 Loss 0.1490 acc 0.7573
        Batchwise validation loss
        Epoch 9 Batch 0 Loss 0.3120 acc 0.7573
        Epoch 9, Loss: 0.24627408385276794, Accuracy: 76.60527801513672, Test Loss: 0.3
        0513450503349304, Test Accuracy: 79.16907501220703
        ===== Start Epoch 10 ======
        Batchwise Train loss
        Epoch 10 Batch 0 Loss 0.2724 acc 0.7570
        Epoch 10 Batch 40 Loss 0.2121 acc 0.7570
        Epoch 10 Batch 80 Loss 0.1736 acc 0.7567
        Batchwise validation loss
        Epoch 10 Batch 0 Loss 0.3214 acc 0.7568
        Epoch 10, Loss: 0.24162065982818604, Accuracy: 76.54540252685547, Test Loss: 0.
        2941944897174835, Test Accuracy: 79.09961700439453
In [ ]: %load_ext tensorboard
In [ ]: | tensorboard --logdir=logs/
        Output hidden; open in https://colab.research.google.com (https://colab.researc
        h.google.com) to view.
        Evaluation
In [ ]: from tensorflow.keras.applications.densenet import preprocess input
        def get_img_tensor(image_path,img_name,model_image):
          img = tf.io.read file(image path+str(img name))
          img = tf.image.decode jpeg(img,channels=3)
          img = tf.image.resize(img,(224,224))
          img = preprocess input(img)
          img_features = model_image(tf.constant(img)[None,:])
          return img_features
```

```
In [ ]: | def evaluate(img name):
          img_tensor = tf.convert_to_tensor([get_img_tensor(img_path,img_name[0], image_f
                                               get_img_tensor(img_path,img_name[1], image_
          img features = tf.constant(img tensor)[None,:]
          feature val = encoder(img features)
          dec_input = tf.expand_dims([tokenizer.word_index['<start>']],1)
          result
          text
          max length output = 60
          for i in range(max_length_output):
            predictions = decoder(dec input, feature val)# the output shape is 1X1X1225
            predictions = tf.reshape(predictions,[predictions.shape[0],predictions.shape[
            # print(predictions)
            # break
            predicted id =tf.argmax(tf.math.log(predictions),1)[0].numpy()
            result.append(tokenizer.index word[predicted id])
            text += " " + tokenizer.index word[predicted id]
            if tokenizer.index word[predicted id] == '<end>':
              return result,text
            dec input = tf.expand dims([predicted id],1)
          return result, text
```

```
In [ ]: import matplotlib.image as mpimg
    from nltk.translate.bleu_score import sentence_bleu
    def test_img_cap(img_data):
        result,text = evaluate(img_data)
        fig,axs = plt.subplots(1,len(img_data),figsize=(10,10),tight_layout=True)
        count = 0
        for img,subplot in zip(img_data,axs.flatten()):
        img_ = mpimg.imread(img_path+img)
        imgplot = axs[count].imshow(img_,cmap='bone')
        count +=1
        plt.show()
        print("Predicted:",text)
```

```
In [ ]: print("Actual",o_test[164])
    test_img_cap(i_test[164])
```

Actual <start> no acute pulmonary disease <end>




Predicted: no acute cardiopulmonary abnormality <end>

In []:

```
In [ ]: print("Actual",o_test[64])
  test_img_cap(i_test[64])
```

Actual <start> no acute cardiopulmonary disease <end>




Predicted: no acute cardiopulmonary abnormality <end>

In []:

```
In [ ]: print("Actual",o_test[66])
    test_img_cap(i_test[66])
```


Actual <start> heart size is normal and lungs are clear stable mm right midlung perform granuloma <end>



Predicted: no acute cardiopulmonary abnormality <end>

```
In [ ]: print("Actual",o_test[266])
   test_img_cap(i_test[266])
```


Actual <start> no acute active cardiac pulmonary pleural disease <end>



Predicted: no acute cardiopulmonary abnormality <end>

```
In [ ]: print("Actual",o_test[29])
   test_img_cap(i_test[29])
```

Actual <start> no acute cardiopulmonary disease <end>



Predicted: no acute cardiopulmonary abnormality <end>

In []:	
---------	--